

Determination of Na_2CO_3 and NaOH in a Mixture by Titration

Aim -

To determine the amount of Na_2CO_3 & NaOH in mixture using HCl .

Principle

When a known volume of mixture of Na_2CO_3 and NaOH is titrated with HCl using Phenolphthalein indicator at its point, all the OH^- ions & only half of the CO_3^{2-} are reacted with HCl acid.

When the titration is continued with methyl orange indicator the remaining half of CO_3^{2-} ions will be neutralised with HCl at the end point.

A = all hydroxide ions of half of carbonate ions

B = half of carbonate ions after phenolphthalein end point.

2B = all carbonate ions

A - B = all hydroxide ions

Titration - I

Sr No.	Vol of Na_2CO_3 (mL)	Burette Reading 1	Burette Reading 2	Volume of HCl (mL)
1	20 mL	0	20.1	20.1
2	20 mL	0	19.5	19.5
3	20 mL	0	19.5	19.5

→ Concordant value = 19.5 mL

Indicator used = Methyl Orange

Calculation - I

Volume of HCl = 19.5 mL

Normality of Na_2CO_3 = 0.05 (N)

Volume of Na_2CO_3 = 20 mL

Normality of HCl = $\frac{0.05 \times 20}{19.5} = 0.0512$ (N).

Procedure -

(i) Titration I: Standardisation of HCl.

Pipette out 20 ml of Na_2CO_3 solution into a clean conical flask and add 2-3 ~~at~~ drops of methyl orange indicator to the soln. Then titrate the solution against HCl acid taken in a burette. Record end point. Repeat titration till concordant value is obtained.

(ii) Titration II: Estimation of Na_2CO_3 and NaOH in a given mixture

Dilute the given unknown solution to 100 ml in a standard flask using diluted water. Pipette out 20 ml of this made up solution into a clean conical flask. Add 2-3 drops of phenolphthalein indicator to the solution & titrate against standardised HCl. Record the burette reading as an end point of the titration when disappearance of pink colour is observed. Let's consider the burette reading at the end point be A ml. To the same soln, add 2-3 drops of methyl orange indicator and continue to titrate till color changes. Note down burette readings as methyl orange

Titration II

Sr No.	Vol of unknown sol	Initial Reading	Burette Readings			
			Vol consumed at Phenolphthalein end point (A)		Vol consumed for methyl orange end point (B)	
			Final	Vol Req	Final	Vol Req
1	20 mL	0	28	28	38.5	10.5
2	20 mL	0	28	28	38.5	10.5

Concordant Values (A) 28 mL (B) 10.5 mL

CALCULATION:

① Estimation of amount of Na_2CO_3

$$\text{Volume of HCl} = 2B = 2 \times 10.5 = 21 \text{ mL}$$

$$\text{Volume of mixture} = 20 \text{ mL}$$

$$\text{Normality of HCl} = 0.05 \text{ N}$$

$$\text{Normality of } \text{Na}_2\text{CO}_3 = \frac{21 \times 0.05}{20} = 0.0525$$

$$\Rightarrow \text{Amount of } \text{Na}_2\text{CO}_3 \text{ present in whole of given solution} = \frac{N_2 \times V_2}{10} = \frac{0.0525 \times 40}{10} = 0.21 \text{ g}$$

② Estimation of amount of NaOH

$$\text{Vol of HCl (V}_1) = A - B = 17.5 \text{ mL}$$

$$\text{Vol of mixture} = 20 \text{ mL}$$

$$\text{Normality of HCl} = 0.05 \text{ N}$$

$$\text{Normality of NaOH (N}_2) = \frac{(A - B) \times N_1}{20} = \frac{17.5 \times 0.05}{20} = 0.04375$$

$$\text{Amount of NaOH present} = \frac{N_2 \times V_2}{10} \left(\text{Eq wt of NaOH} \right) = \frac{0.04375 \times 40}{10} = 0.175 \text{ g}$$

end point. Consider it to be 'B' mL. Repeat titration till concordant value for 'A' and 'B' is obtained.

RESULT -

(i) Amount of Na_2CO_3 present in given solution = 0.2809 g

(ii) Amount of NaOH present in given solution = 0.1792 g